


Folding Suitcase Ramps



Fitting and Maintenance Instructions

CAUTION!

NB: Please read this manual before operation for maintenance instructions and safe usage of these ramps.

Always risk assess any situation for any possible risks.

Only use these ramps with a qualified helper, never exceed the recommended weight limit during use.

Always back the wheelchair/scooter down the ramps and use a lap strap for safety.

Always check the wheelchairs/scooters owners' manual regarding the safe incline for any particular chair.

NB. This equipment must be installed by a competent person.

4ft Folding Ramp

6ft Folding Ramp

Weight Capacity 272 kg (43 st) using both ramps (136 kg (21 stone) per ramp)

Do not exceed the weight limit stated – doing so could put the user at risk.

USING YOUR RAMP

- Always fully extend the ramps before use to minimize the risk of injury, and also make any ascent or descent using the ramps easier.
- Select a place to use where the ramps have the minimum gradient possible. Ensure the gradient of the ramps is not beyond the ability of the wheelchair or occupant.
- Ensure the ramps are safe before use.

TO OPEN & JOIN

- Place the ramp on the narrow edge on the ground, with the handle at the top.
- Pull the locking pins out of the corresponding holes to release the folded ramp.
- Carefully open out the ramp using the hinged join until it forms one continuous ramp.
- Repeat this for the other ramp.
- Position both ramps so that the hinges on the long edges are correctly aligned, then insert the locking pins to fasten the ramps together.

TO CLOSE

- Remove the locking pins from the long edge hinges.
- Carefully fold each individual ramp along the short edge hinge and insert the locking pins into the appropriate hole to fasten the ramp closed.
- Repeat this for the other ramp.

MAINTENANCE & SAFETY

- Regularly check your ramps for worn, loose or damaged parts – If you find any of these to be damaged you must not use your ramps and you must contact your dealer immediately.
- Always use both ramps, and ensure they are level and parallel to each other.
- Ensure ramps are fully extended and locked in position.
- Do not attempt to use ramps on gradients above 1 in 12. Please also check the gradient your wheelchair or scooter is capable of.
- Do not exceed the weight capacity of the ramps.
- Only allow one person/ mobility vehicle to use the ramps at once.
- The surface of the ramps may become slippery in icy or wet conditions. Please exercise caution in these conditions.
- Keep your ramps clean of any dirt or debris as this could interfere with its sliding motion.
- Always use the ramps in line with the wheelchairs/scooter manufacturer's guidelines; NEVER use it at a steeper gradient than recommended.
- Do not modify or add to your ramp as it could void any warranty claims.
- Do not use the ramp for anything other than its original intended purpose.
- Only clean your ramps with mild detergent.
- Use caution when adjusting the length of your ramps and avoid any pinch points.

CLEANING

Clean your Ramps using a mild detergent with a slightly damp, soft cloth only. Only use a mild detergent or disinfectant; never use an abrasive cleaner.

After cleaning, leave the Ramps somewhere warm and dry to air or wipe over with a dry, soft cloth. Do not leave in direct contact with a heat source.