

Shower Chair - Extra Wide



Features

- Wider seat widths and reinforced frame to suit the bariatric client
- Padded seat and curved backrest provide comfort and postural support
- Non-slip rubber feet
- Powder coated steel frame
- Zinc treated to prevent rusting

Specifications

| | |
|-----------------|-------------|
| Seat Height | 470 - 640mm |
| Seat Depth | 450mm |
| Armrest Height | 220mm |
| Back Height | 440mm |
| Max User Weight | 300kg |

Models Available

| | |
|--------|------------------|
| AG0110 | 550mm Seat Width |
| AG0120 | 600mm Seat Width |
| AG0130 | 650mm Seat Width |

