

KING COBRA

Scooter

King
COBRA

*“Powerful performance
and dynamic design”*



CE

Following the success of our Cobra scooter we've developed the new King Cobra. This offers an evolved dynamic design, generous spaciousness, a 31 mile range and a maximum weight limit of 32 stone.

driveTM
www.drivemedical.co.uk

The King Cobra is a revolutionary scooter in terms of design, aesthetics and styling. It has been aerodynamically designed to give an outstandingly stable drive and maximise battery consumption. It's a remarkable scooter that makes a huge statement and sets a new benchmark with its striking dynamic design.



Superb Comfort, Style & Performance...

- » Top speed of 8mph*.
 - » Maximum range of 32 miles on full battery charge*.
 - » Stylish & distinctive design.
 - » Comfortable rotating seat with sliding and reclining facility.
 - » Height adjustable headrest and angle adjustable armrests to suit user requirements.
 - » Electromagnetic brake system and emergency cable brakes.
 - » All round suspension ensures a smooth ride and optimal handling.
 - » Generous tiller adjustment to suit user requirements and assist easy transportation.
 - » User friendly controls with a modern digital LCD display.
 - » Automatic speed reduction when steering for added safety.
 - » Twin USB power socket.
 - » Fully functioning indicators and twin adjustable rear view mirrors for extra visibility.
 - » High level indicators for added safety.
 - » Fitted with reflectors and LED front, rear and brake lights to improve outdoor visibility at dusk or when dark.
 - » Stylish 'Blade' alloy wheels with black pneumatic tyres for a smoother ride.
 - » Delta bars as standard.
 - » Bottle/cup holder included.
 - » Rear anti tip wheels as standard.
 - » Anti roll back system safety device prevents scooter from rolling backwards on hills.
 - » Freewheel facility allows movement of scooter without switching the motor on.
 - » Optional rear pannier box.
 - » Recommended MK Gel batteries: 2 x 12V 96Ah.
 - » Weight capacity 200kg (32st).
- * Subject to conditions of use.



Digital LED display



All round suspension

CODE: **KCOBRAGG**
DESCRIPTION: **King Cobra**



Length:
170cm (67")



Height:
140cm (54")



Width:
73cm (29")



Max Capacity:
200kg (32st)



Ground Clearance:
11.5cm (4.5")



Max Gradient:
10°



Stylish 'Blade' Alloy Wheels



LED Lighting

Leaflet Code: LL208

If you would like to view this product please contact your local stockist, details below.