

PYTHON

Scooter Range

“The Python brings a stylish design and offers a huge 48km range”



The Python range of Scooters are excellent additions to the Drive range, with the performance usually associated with much larger models. These scooters enjoy a whole host of excellent features including all-round suspension, LED headlights and longer range batteries.



www.drive-medical.com.au

drive

DeVilbiss
HEALTHCARE

PYTHON

Scooter Range

Length 121cm (47.5")	Height 122cm (48")
Width 60cm (24")	Weight 100kg (220lb)
Max Capacity 160kg (25st)	Max Gradient 8°



Length 121cm (47.5")	Height 102cm (40")
Width 60cm (24")	Weight 94kg (207lb)
Max Capacity 160kg (25st)	Max Gradient 8°



Length 121cm (47.5")	Height 102cm (40")
Width 60cm (24")	Weight 94kg (207lb))
Max Capacity 160kg (25st)	Max Gradient 8°



PYTHON 8

COLOUR OPTIONS: MS058PBAU-1 – Blue
MS058RDAU-1 –Red

- Max range of 48km*
- Top speed of 13kmph*
- Fully functioning indicators and mirror for road use
- Powerful 500W motor
- Secondary cable brake
- Sport style seat
- Modern dashboard design
- Black pneumatic tyres with large 5" rims
- Hydraulic tiller adjustment
- Delta bars as standard
- Recommended MK batteries: 2 x 12v 40ah**

* Subject to conditions of use. ** Batteries not included.

PYTHON 6

COLOUR OPTIONS: MS051PBAU-1 – Blue
MS051RDAU-1 – Red

- Max range of 48km*
- Top speed of 9kmph*
- Fully functioning indicators and mirror for road use
- All round suspension
- Powerful 470W motor
- Stylish silver wheel hubs
- Delta bars as standard
- Black pneumatic tyres with large 5" rims
- LED front headlight
- Recommended MK batteries: 2 x 12v 40ah**

* Subject to conditions of use. ** Batteries not included.

PYTHON 4

COLOUR OPTIONS: MS050PBAU-1 – Blue
MS050RDAU-1 – Red

- Max range of 48km*
- Top speed of 6kmph
- Powerful 350W motor
- All round suspension
- Stylish silver wheel hubs
- Delta bars as standard
- Black pneumatic tyres with large 5" rims
- LED front headlight
- Recommended MK batteries: 2 x 12v 40ah**

* Subject to conditions of use. ** Batteries not included.