

X² Deluxe Pressure Care Chair



Care Quip introduces the new X² Deluxe Pressure Care Chair. Offering a superior level of pressure care, this is an ideal seating solution for clients sitting for long periods requiring comfort and postural stability. Contoured backrest features lumbar and lateral support and memory foam in seating area contours to the client reducing pressure. 2-way stretch vinyl is antibacterial, vapour-permeable and fluid resistant. The durable and reinforced frame is ideal for facilities and hospitals and features tilt mobility wheels for manoeuvrability.

Features

- High back orthopaedic chair
- Adjustable height and seat depth
- Contoured backrest provides lumbar and lateral support
- Wide padded armrests
- Coffee vinyl upholstery
- 2-way stretch vinyl on seating surface
- Rear tilt mobility wheels
- Powder coated reinforced steel frame

Specifications

- Seat Width: 500mm
- Seat Depth: 460-550mm
- Seat Height: 440-600mm
- Arm Height: 230mm
- Back Height: 750mm
- Overall Width: 700mm
- Max User Weight: 200kg
- Product Code: ED3850

